



# AI. software solutions

NIS-Elements imaging software module

Company names and product names appearing in this brochure are their registered trademarks or trademarks.  
N.B. Export of the products\* in this brochure is controlled under the Japanese Foreign Exchange and Foreign Trade Law.  
Appropriate export procedure shall be required in case of export from Japan.  
\*Products: Hardware and its technical information (including software)

**WARNING** TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. March 2020 ©2020 NIKON CORPORATION



## NIKON CORPORATION

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan  
phone: +81-3-6433-3705 fax: +81-3-6433-3785  
<https://www.healthcare.nikon.com/>

ISO 14001 Certified  
for NIKON CORPORATION

## NIKON INSTRUMENTS INC.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A.  
phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only)  
fax: +1-631-547-0306

<https://www.microscope.healthcare.nikon.com/>

## NIKON INSTRUMENTS EUROPE B.V.

Tripolis 100, Burgenweeshuispad 101, 1076 ER Amsterdam, The Netherlands  
phone: +31-20-7099-000 fax: +31-20-7099-298

[https://www.microscope.healthcare.nikon.com/en\\_EU/](https://www.microscope.healthcare.nikon.com/en_EU/)

## NIKON INSTRUMENTS (SHANGHAI) CO., LTD.

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060  
(Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026  
(Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580

## NIKON CANADA INC.

CANADA phone: +1-905-602-9676 fax: +1-905-602-9953

## NIKON FRANCE S.A.S.

FRANCE phone: +33-1-4516-45-16 fax: +33-1-4516-45-55

## NIKON GMBH

GERMANY phone: +49-211-941-42-20 fax: +49-211-941-43-22

## NIKON INSTRUMENTS S.p.A.

ITALY phone: +39-55-300-96-01 fax: +39-55-30-09-93

## NIKON GMBH SWITZERLAND

SWITZERLAND phone: +41-43-277-28-67 fax: +41-43-277-28-61

## NIKON UK LTD.

UNITED KINGDOM phone: +44-208-247-1717 fax: +44-208-541-4584

## NIKON CEE GMBH

AUSTRIA phone: +43-1-972-6111 fax: +43-1-972-611-140



# AI. software solutions


# Denoise.ai

## Remove shot noise from resonant confocal images


High-speed resonant scanning enables ultra-short dwell times, reducing phototoxicity even during extended periods of timelapse imaging.

However, resonant confocal images contain shot noise as a result of sampling.

With Denoise.ai, the network has been trained to recognize common patterns of shot noise and remove them from resonant confocal data sets. This trained AI algorithm can be used for noise removal in real time. It can greatly improve signal-to-noise ratios (SNR) without requiring line averaging, even when scanning with very short exposure times.



Original




With Denoise.ai

# Clarify.ai


## Remove out-of-focus blur from fluorescence images

During widefield imaging, especially of thick 3D samples, image sharpness may be impaired by blurring caused by signals from sources other than the focal plane.


With Dehaze.ai, the trained AI algorithm allows the user to remove the out-of-focus blur from the image in real time. It enables acquisition of clear, haze-free images from deep within a sample without requiring traditional deconvolution even when signals are weak, and greatly reduces photo toxicity.



Original



With Clarify.ai



- Provides bright, uniform illumination across an unprecedented 25 mm field of view that maximizes the sensor area of large-format CMOS cameras and significantly improves data throughput

- The hardware-triggering capabilities of Ti2-E enhance even the most challenging, high-speed imaging applications

Motorized Inverted Microscope

## ECLIPSE Ti2-E

Confocal Microscope

## A1R HD25


- When used with the Ti2-E, the A1R HD25 is capable of 25mm FOV images, enabling increased screening analysis throughput and large image stitching
- Equipped with both a resonant scanner capable of a maximum ultrafast imaging speed of 720 fps and a galvano scanner capable of a maximum imaging resolution of 4096 x 4096 pixels.



# Enhance.ai

## Restores details in under-exposed images


Imaging of dim fluorescent samples or applications that require low-level light exposure typically result in compromised image quality with poor signal-to-noise ratio. Enhance.ai learns what a high signal-to-noise image looks like by comparing under-exposed and optimally-exposed images. Enhance.ai can then restore details in under-exposed or dim fluorescent images, enabling researchers to gain more insights from their low signal imaging applications.




Original

Enhance.ai

DAPI stained nuclei are purposefully underexposed to limit the specimen's exposure to near-UV light. Enhance.ai is used to restore the signal-to-noise ratio to normally exposed DAPI staining, for easy segmentation and counting.



Original



With Convert.ai

DAPI staining of nuclei is a common method allowing cell counting and segmentation. Convert.ai can be trained to predict where the DAPI label is present in DIC or phase images. This predicted channel can then be used for segmentation and counting, without ever having to label the specimen with DAPI or acquire a fluorescence channel.

# Convert.ai


## Predicts DAPI label locations based on unstained images

Convert .ai learns associated patterns present in two different imaging channels. After training, Convert.ai can then predict the pattern in the second channel even when presented with only the first channel. DAPI-based staining of nuclei is a common method for cell segmentation and counting. Convert.ai can be trained to predict where the DAPI label would be based on unstained DIC or phase-contrast images. This enables users to perform nuclei-based image analysis without having to stain samples with DAPI or acquire a fluorescent channel.


# Segment.ai

## Segments structures that cannot be identified by traditional thresholding

Segment.ai enables complex structures to be easily identified and segmented. Neurites in phase-contrast images are traditionally difficult to define by classic thresholding. Segment.ai can be trained with a small subset of hand-traced neurites to automatically detect and segment neurites from thousands of untraced datasets.



Original



Segment.ai

Neurites in phase-contrast were not possible to define accurately by traditional thresholding. Segment.ai was trained on hand-traced neurites (human recognized) and learned how to trace neurites in subsequent images.

## Expands the NIS-Elements platform by incorporating tailor-made solutions for acquisition, visualization and analysis

Enhance.ai/Convert.ai/Segment.ai employs convolutional neural networks (CNNs) to learn from labeled training data created by either conventional segmentation or human-assisted tracing of small subsets of representative samples. When using the module, the software interface makes it easy to apply complex deep learning to sample data, eliminating the need to design a complex neural network and apply training data to it.

Automated tools take this training data and apply it to the neural network to recognize patterns. The resulting training recipe can then be applied repeatedly and reliably to similar samples to process or analyze huge volumes of data at significantly faster speeds than traditional techniques.

